
en	Instructions for use	hr	Upute za upotrebu
da	Brugsanvisning	cs	Návod k použití
nl	Instructie voor gebruik	hu	Használati útmutató
fi	Käyttöohjeet	no	Bruksanvisning
fr	Mode d'emploi	pl	Instrukcja użytkowania
de	Gebrauchsanleitung	ro	Instrucțiuni pentru utilizare
el	Οδηγίες χρήσης	sk	Návod na použitie
it	Istruzioni per l'uso	tr	Kullanma Talimatı
pt	Instruções de Uso	zh	使用说明
ru	Инструкции по применению	ja	使用方法
es	Instrucciones de uso	ko	사용 지침
sv	Bruksanvisning till användning	et	Kasutamisujuhend
bg	Инструкции за употреба	lt	Naudojimo instrukcija

CritiCath™ Flow-Directed Thermo dilution Catheters Models - SP5105, SP5107, SP5127, SP5507, SP5527 Instructions

Models with Suffix :
 "H" are heparin complex coated.
 "L" are longer (110 cm) in SP5105.
 "S" are softer, more flexible catheters.
 "U" are polyurethane (body only).
 "P" are polyurethane (entire catheter except balloon).
 "13", "14", "23" and "24" have pre-attached protection sheaths; "23" & "24" have SP5045 also.
 "M" have port locations for smaller anatomies.

CONTRAINDICATIONS

Use of a heparin complex coated catheter is contraindicated in patients with known sensitivity to heparin or benzalkonium chloride.

General Description

The CritiCath™ Flow-Directed Catheters are 14kΩ thermomodulation catheters designed for use as a diagnostic tool. Model SP5105 is designed for use in pediatric patients while Models SP5105L, SP5107, SP5127, SP5507 and SP5527 are for adult patients. These catheters allow physicians to measure intracardiac pressures and cardiac output when used with Argon Critical Care Systems (Formerly known as BD Critical Care Systems) or other similar Cardiac Output Computers and appropriate pressure monitoring equipment. The catheters are intended for use at bedside, in cardiac catheterization laboratory, surgery suite, post anesthesia recovery unit, and other specialized critical care units. In addition to measuring cardiac output (thermodilution method) these catheters may be used to measure right atrial (RAP), pulmonary artery (PAP) and pulmonary capillary wedge pressures (PCWP), to continuously monitor PA temperature, to withdraw blood samples from the RA and PA, and to administer drugs and I.V. solutions.

The catheter body material is Polyvinyl Chloride unless indicated by U or P which indicates polyurethane.

These catheters are designed for diagnostic procedures and do not require the use of fluoroscopy for insertion. However, they are radiopaque so that fluoroscopy can be used to guide insertion and to verify position following insertion. The catheters are marked every 10 cm to indicate distance from the distal tip. Narrow bands represent 10 cm markings and wide bands represent 50 cm markings.

SP5105, SP5107 and SP5127 has 4 lumens :

- **'PA DISTAL' [Distal (PA) Lumen]:** Terminates at catheter tip and is used for measurement of PA, PCWP, mixed venous pressure and for blood sampling.
- **'CVP PROXIMAL' [Proximal (CVP) Lumen]:** Terminates at 15 cm (for SP5105), 29 cm (for SP5105L, SP5107, SP5507, SP5127, SP5527) and 25 cm (for "M" version models) proximal to catheter tip and is used for injecting thermal bolus into the RA, for measuring RAP or CVP and for blood sampling.
- **'BALLOON .75 mL CAP'/'BALLOON 1.5 mL CAP' [Balloon Lumen]:** Terminates in 0.75mL balloon (for SP5105 models) and 1.5mL balloon (for the other models) near catheter tip and allows inflation and deflation of the balloon during catheter insertion and measurement of PCWP.

- **Thermistor Lumen:** Terminates in a thermistor bead 1.5 cm and 3.5 cm proximal to the tip for pediatric & adult respectively. Allows cardiac output computation and PA blood temperature measurement.

SP5127 is designed for femoral insertion (Use in cardiac catheterization laboratory.)

SP5507 and SP5527 contains five lumens. Lumen design is as indicated above with an additional lumen:

- **'CVP MEDICATION' [Medication (CVP) Lumen]:** Terminates 31cm (27cm for "M" version catheters) from catheter tip. May be used for solution infusion, cardiac output injections, blood sampling, fluid and drug administration and monitoring of RAP.
- **Ventricular Lumen:** Terminate 17cm from catheter tip. Used mainly for monitoring of right ventricular pressure (RVP) for SP5527 only.

SP5527 is a ventricular catheter.

NOTE: Administration of blood through the 5-lumen catheter is not recommended.

A damaged catheter cannot be repaired. The packages are designed to prevent catheter kinking and to protect the fragile balloon and heparin complex coated surface.

Catheter Testing

1. Use aseptic method to remove catheter.
2. (Models "13", "14", "23", and "24" Only) Tighten the Tuohy-Borst adaptor on the protection sheath securely to the catheter.
- NOTE:** Do not overtighten the adaptor as this may occlude catheter lumens and result in inaccurate pressure readings.
3. Test thermistor by connecting catheter to Cardiac Output Computer. If computer's display indicates appropriate room temperature, the thermistor is operating properly.

CAUTION: Do not use a standard continuity checker or ohmmeter when the catheter is in the vascular system. The relatively high current in the meter can cause electric shock to the patient in the event of insulation breakdown. Excessive manipulation of catheter tip may damage thermistor and result in intermittent or no output readings.

4. Test balloon by placing it in sterile water and injecting 0.75mL of air for SP5105 and SP5105L (1.5mL for SP5107, SP5127, SP5507, SP5527). If air bubbles appear around balloon, do not use the catheter. Deflate balloon by removing syringe.

CAUTION

- Never use liquid for balloon inflation.
- Always deflate the balloon by removing syringe; balloon may be damaged by aspirating air with syringe.

NOTE: Liquids within balloon lumen may cause balloon to remain inflated. Special care must be taken to prevent moisture from entering balloon lumen.

5. Evacuate air from injectate and IV solution bags.

NOTE: Eliminating air from solution will prevent air from entering the system when solution is exhausted or bag inverted.

6. Connect stopcocks to the ends of the CVP proximal and PA distal lumens. Flush and fill the CVP proximal and PA lumens with sterile solution.

WARNING: Stopcock lever must be positioned 90° for "OFF" position. Do not position stopcock at 45° to achieve an "OFF" position. 45° positioning of stopcock is imprecise and may result in contamination, bleed back, or air embolism.

Catheter Insertion

Catheter Insertion should follow the basic procedure recommended below. However, catheter use must always conform to your hospital policies and procedures. Thorough understanding of listed "Cautions" is helpful before use of these catheters.

1. Monitor ECG continuously during catheter insertion.
2. Test catheter as outlined in "Catheter Testing." Connect stopcocks to the CVP and PA lumens. When using an In-Line Injectate Temperature Sensor, connect it directly to CVP lumen, and then add the stopcock. Fill both lumens with sterile solution, to remove air bubbles. Connect PA distal lumen to a pressure-monitoring system.
3. Insert catheter percutaneously or via a cutdown. Typical insertion sites include the median basilic, jugular, femoral and subclavian veins.
4. Advance catheter into vena cava and partially inflate the balloon to 0.5 mL (for SP5105 & SP5105L) and 1 mL (for SP5107, SP5127, SP5507 & SP5527). Filtered CO₂ is recommended for inflation because of its rapid absorption into blood in case of balloon rupture. However, room air is frequently used if there is no right-to-left shunt or pulmonary arteriovenous fistula. Introduction of an air embolus into the arterial system can cause serious complications. The risk of balloon rupture and likelihood of air entering the arterial system must be considered when selecting air as the inflation medium. Determine catheter length inserted by referring to 10 cm interval markings on catheter. Follow the pressure waveforms as the catheter is advanced (reference Figure 1).
5. When RAP waveform is seen (Figure 1-A), inflate balloon to full capacity (0.75 mL for SP5105 & SP5105L and 1.5 mL for SP5107, SP5127, SP5507 & SP5527).
6. Advance catheter through RA and into RV (Figure 1-B). If RVP tracing is not recorded after the catheter has been advanced beyond RA, deflate balloon by removing syringe, pull catheter back slowly, reinflate and advance catheter (Figure 1-C). Observe a RVP waveform. Advance into PA. If the catheter does not pass into PA, deflate balloon and withdraw the catheter until tip is in RA. Then reinflate balloon and repeat steps 4, 5 and 6.


Figure 1A-E Pressure Waveforms During Catheter Insertion and Positioning. (Chart Speed: 25 mm/sec.)

NOTE: If balloon has been inflated for more than 4 or 5 minutes during insertion, the inflated balloon volume may decrease due to air diffusion through latex balloon. Remove syringe and reinflate balloon if necessary.

7. Continue advancing the catheter slowly through PA until a PCWP is seen. (Figure 1-D).

8. Check proper positioning of catheter in PA by deflating balloon completely (removing syringe) and observing the change from a PCWP tracing to PAP tracing (Figure 1-E). Gradually, reinflate balloon until a PCWP is obtained and record the volume of air required. It should take 0.5 to 0.75mL (for SP5105 & SP5105L) and 1.25 to 1.5mL (for SP5107, SP5127, SP5507, SP5527). If less than 0.5mL (for SP5105 & SP5105) or 1.25mL (for SP5107, SP5127, SP5507 & SP5527) is required, the catheter may be advanced too far, increasing the likelihood of distal migration and spontaneous wedging. Pull catheter back 2 or 3 cm and recheck inflation volume.

9. (Models "13", "14", "23", and "24" Only) Extend the protection sheath over catheter body and fasten to hub of introducer.

Cardiac Output Measurements

For cardiac output measurement instructions, refer to operating manual provided with your cardiac output computer. Suggested introducer sizes for catheters :

Catheter	Introducer / With Sideport Flow	
5F	6F	
7F	7.5F	8F
7.5F	8F	8.5F

Model SP5105 is specifically designed to operate with the following computation constants :

	Volume Injected mL	Computation Constant
• Indicator at Ice	5	.291
Temperature, In-Line	3	.170
• Indicator at Ice Bath	5	.279
Temperature (0-5°C)	3	.160
• Indicator at Room	5	.316
Temperature (23 to 25°C)	3	.188
(the same for 18 to 22°C)		

Models SP5107, SP5127, SP5507 and SP5527 uses the following computation constants:

	Volume Injected mL	Computation Constant
• Indicator at Ice	10	.579
Temperature, In-Line	5	.281
	3	.160
• Indicator at Ice	10	.566
Bath Temperature (0-5°C)	5	.270
	3	.151
• Indicator at Room	10	.628
Temperature	5	.309
(23 to 25°C)	3	.181
(the same for 18 to 22°C)		

Cautions

1. To avoid damage to catheter or balloon when inserting the catheter by the cutdown procedure, it is recommended that a vessel dilator or disposable vein guide introducer be used. **NEVER** use forceps on the catheter.
2. (For Heparin Complex Coated Catheters only) If catheter surface is treated with heparin complex coating, do not wipe the catheter. The percutaneous method of insertion can minimize the possibility of bleeding at the entry site. If wound closure is prolonged, an appropriate topical hemostatic agent should be used to neutralize the effect of any heparin complex that may have wiped off at the edge of the wound.
3. It is generally recommended that the catheter not be left in the patient for longer than 72 hours.
4. Catheters should not be advanced after properly positioned as the portion of catheter left outside the body may not be sterile. If a catheter protective sheath is used, catheter may be repositioned as needed.
5. Never make a cardiac output measurement with the balloon inflated.
6. Determine PCWP by inflating balloon slowly, stop when PAP waveform changes to PCWP waveform. Deflate balloon after completing measurement.
7. To minimize ventricular irritation, always inflate the balloon before the catheter reaches RV.
8. Do not exceed the recommended 0.75mL (for SP5105 & SP5105L) and 1.5mL (for SP5107, SP5127, SP5507 & SP5527) balloon volume.
9. Initial placement of catheter in PA should always be made with the maximum recommended balloon inflation volume. An under inflated balloon will be smaller, thus allowing catheter to be positioned in a narrower portion of PA, increasing the likelihood of spontaneous wedging or the risk of PA rupture.
10. A slow heparin drip may be used to maintain catheter patency and should be implemented before catheter insertion and maintained thereafter. Be sure to backflush lumen after drawing blood. Care must be taken not to damage the balloon when advancing catheter through hemostasis valves, catheter/sheath adapters and introducers. Best results are obtained by dilating valves, rubber seals and moistening surfaces of balloon seals and valves with a sterile, injectable fluid. A flow-directed catheter may migrate into the distal PA and cause spontaneous wedging. To detect wedging, PA pressure waveform should be monitored continuously or at short intervals.

Complications

It is important for physicians to be familiar with the catheter insertion procedures, its advantages and associated risks in order to reduce the incidence of complications. Risks associated with the use of this device include: Dysrhythmias, Pneumothorax, Cardiac Tamponade, Balloon Rupture, Thrombus/ Pulmonary Infarct, Perforation of the Pulmonary Artery, Knotting, Sepsis/Infection, and other complications including right bundle branch block (RBBB)

For further information regarding complications, contact your local representative.

**REFERENCE ONLY
DOCUMENT CONTROL**

Storage

Store catheter in a cool, dark area away from fluorescent lights and sunlight to prevent premature deterioration of the latex balloon.

STERILE and non-pyrogenic in unopened, undamaged package. For single use only. Check integrity of the individual package before use. Dispose of product after use. Do not resterilize.

Re-use may lead to infection or other illness/injury.

This product contains a phthalate which may cause a risk to (children (male), pregnant or nursing women) in conjunction with prolonged duration of use.

For reordering information or assistance please contact local representative.

Critiath™ Flow-Directed Thermodilution Catheters Models - SP5105, SP5107, SP5127, SP5507, SP5527

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Argon Critical Care Systems Singapore Pte. Ltd.,
Argon Medical Devices Inc., 1445, Flat Creek Road, Athens,